

Reference class systems for experts in the field of  
**PSYCHIATRY & PHYSIOLOGY**

**RMS QUADRAMIND**

(Multibehavior Therapy)



RMS QUADRAMIND (Multibehavior Therapy Unit) is Combination of Sex Therapy, Electrosleep, Brain Polariser and Aversion Therapy in a single unit.



**RECORDERS & MEDICARE SYSTEMS (P) LTD.**

**Certifications: ISO 9001:2008, EN ISO 13485:2003**



All the four parameters work independently and can be used Simultaneously.

1. **Sex Therapy:-** Sex Therapy instrument is used for the treatment of sexual dysfunctions arising due to psychological reasons.

**Features & Specifications**

- Male & Female therapy can be used Simultaneously.
- Stimulus Intensity : Variable for Male 0 to 10 volts.
- Vibrating intensity : Variable for Female 0 to 5 volts.
- Stimulating Electrode : Male - Rubber pad electrode
- Vibrating electrode : Female - Vaginal shaped electrode

2. **Electrosleep:-** It is commonly used in the treatment of chronic insomnia, anxiety and depression but has also been used for psychosomatic disorders such as asthma, tension headache and for organic disorders including hypertension.

**Features & Specification**

- Consists of the application of short duration, low amplitude pulses of direct current to the patient's brain via externally placed electrodes.
- Causes natural physiological sleep normalizing activity of nervous systems, psychoemotional activity and functioning of internal organs.
- Stimulus Intensity :- 0 to 30 v Variable
- Audio Stimulus :- Gives audio stimulus along with electrical impulses.
- Based on SOMLEC therapy.

**Electrodes:-** Velcro tape electrodes.

3. **Brain Polariser:-** The Brain Polariser instrument is used for the treatment of depression . It is a fully non convulsive therapy. It is beneficial to the patients even where ECT & drugs fail.

**Features & Specification**

- Treatment is provided with the pulsating low intensity current.
- Stimulus Intensity :- 0 to 10V variable.

**Electrodes :-** Velcro tape electrodes.

4. **Aversion Therapy:-** This instrument is used to administer low current electrical shock or other naxious stimuli to the patient to modify undesirable behaviour characteristics like addicts from alcoholism, chain smoking etc.

**Features & Specifications**

- Stimulus Intensity - 0 to 100V variable
- Operation through remote control for stimulation.

**Electrodes:-** Velcro tape electrodes.

Since R&D is a continuous process features & specifications are subject to change without notice.



**Recorders & Medicare Systems (P) Ltd.**

**H.O. & Works:** 196, Industrial Area, Phase 1, Panchkula-134113 (Haryana), INDIA.

Ph.: +91-172-2564196, 2565196, Fax: +91-172-2566196

website: [www.rmsindia.com](http://www.rmsindia.com), e-mail: [helpdesk@rmsindia.com](mailto:helpdesk@rmsindia.com)

**Regd. Office:** 181/5, Industrial Area, Phase 1, Chandigarh-160002, INDIA